



# **Blood Chemistry Sensors**

Suresh Devasahayam  
Department of Bioengineering  
Christian Medical College, Vellore

# Outline

---


- pH sensors
- Oxygen sensors
- Glucose sensors
- CO<sub>2</sub> sensors

# Normal Ranges


- $\text{PO}_2$ : 80-104 mmHg
- $\text{PCO}_2$ : 33-48 mmHg
- pH : 7.31-7.45
- Glucose: 70-110

# Photospectrometry

- Estimate the spectrum of optical properties of a sample
- Measure intensity of light at few selected wavelengths



# Optical Spectrum estimation



$$n \lambda = d \sin \theta$$

$d$  = grating spacing


# LED colours

## LED Color Guide


LED P/N Suffix	Description	Chemistry	# of Elements	Color Temperature (CCT Typ)	Peak Wavelength (λ / x-coord)	Dominant Wavelength (λ / y-coord)	Forward Voltage (Vf Typ) (Vf Max)		Brightness
H	High Efficiency Red	GaP	2	~	700	660	2.0	2.5	Standard
SR	Super Red	GaAlAs	3	~	660	640	1.7	2.2	High
SR	Super Red	AlInGaP	4	~	660	640	2.1	2.5	High
SI	Super High Intensity Red	AlInGaP	4	~	636	628	2.0	2.6	High
I	High Intensity Red	GaAsP	3	~	635	625	2.0	2.5	Standard
ZI	TS AlInGaP Red	AlInGaP	4	~	640	630	2.2	2.8	High
SO	Super Orange	AlInGaP	4	~	610	602	2.0	2.5	Standard
A	Amber	GaAsP	3	~	605	610	2.0	2.5	Standard
SY	Super Yellow	AlInGaP	4	~	590	588	2.0	2.5	Standard
ZY	TS AlInGaP Yellow	AlInGaP	4	~	590	589	2.3	2.8	High
Y	Yellow	GaAsP	3	~	590	588	2.1	2.5	Standard
SUG	Super Ultra Green	AlInGaP	4	~	574	568	2.2	2.6	High
G	Green	GaP	2	~	565	568	2.2	2.6	Standard
SG	Super Green	GaP	2	~	565	568	2.2	2.6	Standard
PG	Pure Green	GaP	2	~	555	555	2.1	2.5	Standard
UPG	Ultra Pure Green	InGaN	3	~	525	520	3.5	4.0	High
UEG	Ultra Emerald Green	InGaN	3	~	500	505	3.5	4.0	High
USB	Ultra Super Blue	InGaN	3	~	470	470	3.5	4.0	High
UV	Ultra Violet	InGaN	3	~	410	~	3.5	4.0	Standard
SUV	Super Violet	InGaN	3	~	380	~	3.4	3.9	Standard
T	Turquoise	InGaN	3	~	0.19	0.41	3.2	4.0	Standard
V	Violet / Purple	InGaN	3	~	0.22	0.11	3.2	4.0	Standard
P	Pink	InGaN	3	~	0.33	0.21	3.2	4.0	Standard
MW (Warm)	Warm White	InGaN	3	3000K	~	~	3.3	4.0	High
NW (Neutral)	Neutral White	InGaN	3	4000K	~	~	3.3	4.0	High
UW (Cool)	Cool White	InGaN	3	6000K	~	~	3.3	4.0	High



# LED spectra




# White light using phosphor coated blue LEDs







# Absorption spectrum change of phenol red with pH




# Peak spectrum shift of phenol red with pH




## $O_2$ fluorescence dye (e.g., Anthracenecarboxaldehyde, $C_{15}H_{10}O$ )




# Light absorption by tissue




# Light absorption and oxygenated haemoglobin



# Ion-sensitive Field Effect Transistors



# ISFET with sample holder





# Intravascular sensor

